



DENTIRO® Sensitive

Surface disinfectant with a pleasant lemon scent suitable for use on almost all surfaces

- Pleasant lemon fragrance
- Excellent material compatibility
- Low alcohol content of less than 50%

DENTIRO® Sensitive is an aldehyde-free, ready-to-use product based on alcohols for the quick disinfection and cleaning of surfaces of non-invasive medical devices (e.g. laboratory benches, treatment chairs, unit surfaces, instrument trays, hand-pieces, contra-angles, operating light handles, and dental instruments made of steel, aluminum, copper, silicon, or rubber). The low alcohol content of less than 50% ensures high material compatibility and reduces the potential for allergies and irritation. DENTIRO® Sensitive is easily biodegradable and has a discrete lemon fragrance. DENTIRO® Sensitive can be used to remoisten DENTIRO® Wipes and ISORAPID® Wipes. To make your own pre-moistened wipes, use DENTIRO® Dry Wipes or ISORAPID® XXL Wipes.

Use



Completely moisten the surface or object to be disinfected and leave to act. Adhere to exposure times. If necessary, dry the surface or object with a clean single-use cloth. Use on alcohol compatible materials only. Test before use on sensitive materials by applying to a small inconspicuous area.

Scope of Application



- For the fast-acting spray and wipe disinfection as well as the cleaning of surfaces of non-invasive medical devices, e.g. laboratory benches, treatment chairs, unit surfaces, instrument trays, hand-pieces, contra-angles, operating light handles, and instruments made of steel, aluminum, copper, silicon, or rubber

Precautions



R10 Flammable. R41 Risk of serious damage to eyes. S2 Keep out of the reach of children. S16 Keep away from sources of ignition – No smoking. S23 Do not breathe spray. S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. S51 Use only in well-ventilated areas. Use disinfectants safely. Always read label and product information before use.

Delivery Forms



100 ml with clip lid, 500 ml with clip lid, 1 liter with clip lid, 5 liter, 10 liter



Effective Spectrum / Exposure Times

Bactericide / Fungicide (<i>C. albicans</i>) (1)	
(clean conditions)	30 seconds
(dirty conditions)	60 seconds
Tuberculocide (<i>M. terrae</i>) (1)	60 seconds
Mycobactericide (<i>M. avium</i>) (1)	60 seconds
Hepatitis B / HIV (2)	30 seconds
BVDV / Hepatitis C (3)	30 seconds
Herpes simplex virus	15 seconds
Influenza A virus (incl. H1N1)	60 seconds
Adenovirus (4)	60 seconds
Polyomavirus (SV 40) (4)	60 seconds

1) DGHM Standard Methods for surface disinfection incl. tests under practical conditions; 2) HBsAg antigen test; 3) BVDV surrogate virus test; 4) Non-enveloped virus

Complies with DGHM Standard Methods and EN 1040, EN 1276, EN 1275, EN 1650, EN 14348; DGHM / VAH listed



Composition in 100 g

20 g Ethanol, 28 g 1-Propanol, 0.056 g Quaternary ammonium compounds.